#### Step 3

We use thin glass spotting tubes to "spot" our TLC plate. It is important to not overload or underload the place. For this reason, we use a dilute solution of analyte and check after each application under a UV lamp.

Lighht

### Step 1

We begin with a blank TLC plate.

Predominantly, these are Silica (SiO<sub>2</sub>), a polar inorganic powder, bound to plastic or glass.

Eluent moves up the TLC plate by capillary action.

As it does, the analytes spotted onto the plate equilibrate between the eluent, or mobile phase,

and the silica, or stationary phase.

Silica (SiO<sub>2</sub>) is exceedingly polar compared to the eluent in every case.
As such, **less polar compounds** will spend

more time moving within the eluent mobile

More polar compounds will spend more time stationary in the stationary silica phase.

Therefore, less polar compounds move further up the plate than more polar compounds do.

Solvent System

# Step 5 You will notice we cannot directly see where the analytes have moved to on the plate. We visualize the TLC plate using UV light; many compounds with pi systems (double bonds, etc) will absorb UV light. Because the silica is fluorescently doped, Step 7 wherever an analyte that can absorb UV light is will appear dark. We use compounds' RF values to make We circle the dark spots with a pencil so that predictions about their polarities and, we may see them once the UV light is if we have a standard to compare unknown removed. analytes to, their identities. Solvent System RF = 0.2RF = 0.0.5 most polar |---|----|----|----|----|----|

Step 6

#### Step 2

chamber.

We mark the origin with a line and place tick marks for each analyte. It is important to place the origin high enough on the plate such that it does not go below the solvent in the

Solvent System

We also write which solvent system we will use near the top of the plate.

This is done with penil. Pen ink would dissolve and separate!


#### Step

We place the loaded TLC plate into the developing chamber containing solvent or *eluent*.

## Once the solvent front nears the top of the plate, we carefully remove the TLC plate from

the developing chamber,

UV light will not be absorbed by all compounds, so we also use a *stain* (in this case, lodine) to chemically alter the spots such that they can be seen with the naked eye colorimetrically.



The analyte equilibrates between the mobile

(solvent) and stationary (silica) phases

Plate is placed into chamber with solvent

The analyte equilibrates between the mobile

(solvent) and stationary (silica) phases

